


TL-2100G Cordless Fall Monitor



Programming/ Setup

- ◆ Remove protective boot, place 3 x AA batteries in rear & adjust volume & tone. Replace boot.
- ◆ Place 1 or 2 sensor pads (switched on) within reach. Press & hold programming button on rear until you hear 3 beeps. Release the button – you are now in programming mode.
- ◆ Pad 1 light on monitor will illuminate, place pressure on 1st pad for 3 seconds then release. Pad is now programmed as pad 1
- ◆ Pad 2 light on monitor will illuminate, place pressure on 2nd pad for 3 seconds then release. Pad is now programmed as pad 2.
- ◆ Press reset button on front of unit to exit programming (If programming just one pad, press reset button twice to exit programming)
- ◆ Testing – apply pressure to pad 1, monitor will beep twice. Status light for pad 1 will blink while pressure remains on pad. Remove pressure. alarm will sound. Press reset button to silence alarm (or use caregiver key)

Notes

- ◆ For use with one OR two cordless bed pads, chair pads, floor mats.
- ◆ 5m maximum distance between monitor & sensors
- ◆ Volume dial, tone settings inside battery case at rear of monitor
- ◆ Reset switch on side of monitor – "Reset Switch On" means the reset button on front of monitor is active. "Reset Switch Off" means Caregiver Key point at front of monitor is active, and reset button on front is inactive.
- ◆ When programming a floor mat, start with the floor mat off, then switch on to programme into monitor (Equivalent to placing pressure on the bed/ chair pads)

