




TL-2100E Economy Fall Monitor



Programming/ Setup

- ◆ Remove protective boot, place batteries in rear & adjust volume & tone. Replace boot.
- ◆ Connect the corded pressure pad (socket located at bottom of monitor – run cable through strain relief recess)
- ◆ Testing – Slide power switch on side of monitor to on. Apply pressure to sensor pad to activate. The monitor will beep twice to indicator pad & monitor are in use. Status light will blink every 4 seconds while pressure remains on pad. When pressure is removed from pad, alert tone will sound. Slide power switch on side of monitor to off to silence/reset. Slide switch to on to reactivate the system

Notes

- ◆ For use with one corded bed pad, chair pad, floor mat or seat belt.
- ◆ Volume dial, tone settings inside battery case at rear of monitor
- ◆ When programming a floor mat, start with the floor mat off, then switch on to programme into monitor (Equivalent to placing pressure on the bed/chair pads)

